



Create a new era of rapid molecular diagnosis

iPonatic

Product features

Fast	Simple	Easy	Intelligence
Get result within 15 minutes	Pre-packaged test strip Room temperature preservation Opening while using	Easy to use Integration of Extraction and Amplification	Intelligent voice prompt Analyse results automatically

Test Results

Expert Method Start Amplification Real-time state **Data** Standard Curve Return

Test list :

Module list : Channel list :

FAM
 VIC
 ROX
 CY5

Amplification curve Melting curve Standard curve

Binding curve

Real time fluorescence signal

Amplification cycle

Data analysis of hundreds of samples shows that linear correlation R2 between iPonatic and manual operation reached 0.9667. At the same time, the systematic bias of the two groups of data is only ±0.01 log10. It means there is good consistency between iPonatic and manual operation.

